

3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

# **Product Information**

# CompoZr® Disease Model Cell Lines MCF10A Cells SYK -/-

Catalog Number **CLLS1048**Storage Temperature –196 °C (liquid nitrogen)

# **Product Description**

CompoZr® zinc finger nuclease (ZFN) technology is a fast and reliable way to manipulate the genome in a targeted fashion. ZFNs are synthetic proteins engineered to bind DNA at a sequence-specific location and create a double strand break (www.compozrzfn.com). The cell's natural machinery repairs the break in one of two ways: non-homologous end joining or homologous recombination. The non-homologous end joining pathway resulted in deletions at the SYK locus (see Figure 1). Single cell knockout clones were isolated and followed for more than twenty passages to establish stable cell lines.

While the targeted gene in this cell line is diploid, ZFN-mediated gene knockout technology is not limited to diploid targets, allowing the researcher to pursue many of the polyploid cell lines often characteristic of cancer. Modified cell lines provide the basis for the development of various assays for compound screening. Here, the target gene and resulting protein are disrupted, in contrast to cell lines with normal expression.

The human gene known as spleen tyrosine kinase (SYK) is located on chromosome 9 and encodes a non-receptor tyrosine kinase. This kinase was first recognized as a 40 kDa proteolytic fragment derived from a p72 tyrosine kinase present in various hematopoietic cells. SYK has been recognized as a tumor supressor where the level of expression has been inversely correlated with malignancy. Expression of SYK is thought to help maintain breast tissue integrity by promoting cell-cell adhesions.

Studies employing RNAi to knockdown expression levels in near-normal breast cell lines have demonstrated loss of SYK leads to loss of vimentin and promotes the epithelial to mesenchymal transition (EMT).<sup>4</sup> SYK has been associated with inhibition PI3K and NF<sub>K</sub>B,<sup>5</sup> as well as playing a role in the phosphorylation of E-cadherin, demonstrating further its significance in cell-cell adhesion.<sup>6</sup>

*In situ* loss of expression of SYK is generally associated with methylation of the promoter, <sup>7</sup> and lowered expression has been associated with an increased risk of developing distant metastasis. <sup>8</sup>

For further information and to download sequence of modified locus, go to the website: www.wherebiobegins.com/biocells

### Figure 1.

Creation of SYK Knockout in MCF10A Cells

Site-specific deletion at the SYK Locus in MCF10A cell line, Alleles 1 and 2 – 22 bp deletion:

Schematic of the genomic sequence at the target region (exon 2) recognized by the ZFN pair; the resulting deletion, and the CEL-I primer sequences:

CEL-I Primers - <u>Bolded and underlined</u>
ZFN binding site - <u>UPPER CASE</u>, <u>BOLDED RED</u>
ZFN cut site - lower case red
Deletion - <u>yellow highlighted</u>

Genotype: del 22/del 22 (homozygous)

#### Components

MCF10A mutant cell line with SYK gene 1 vial

knocked out

Catalog No. CLL1048

Parental mammary epithelial cell line 1 vial (ATCC<sup>®</sup> Catalog No. CRL-10317<sup>™</sup>)

Catalog No. CLL1040

#### Cell Line Description

1 vial of modified MCF10a cells contains  $\sim$ 2 × 10<sup>6</sup> cells.

Organism: Homo sapiens (human)

Tissue: mammary gland; breast

Age: 36 years

Gender: Female

Ethnicity: Caucasian

Morphology: Epithelial

Growth properties: Adherent

DNA profile

Short Tandem Repeat (STR) analysis:

Amelogenin: X CSF1PO: 10,12 D13S317: 8,9 D16S539: 11,12 D5S818: 10,13 D7S820: 10,11 THO1: 8,9.3 TPOX: 9,11 vWA: 15,17

Parental Cell Line: ATCC Catalog No. CRL-10317

Note: Please see CRL-10317 product datasheet from ATCC for additional information about the origin of these cell lines. Cytogenetic information is based on initial seed stock at Sigma Life Science. Cytogenetic instability has been reported in the literature for some

cell lines.

Complete Medium: Dulbecco's Modified Eagle's Medium (DMEM)/Ham's Nutrient Mixture F12 (1:1) with 2.5 mM L-glutamine, 5% horse serum, 10  $\mu$ g/mL human insulin, 0.5  $\mu$ g/mL hydrocortisone, 10 ng/mL EGF, and 100 ng/mL cholera toxin. This medium is formulated for use with a 5% CO<sub>2</sub> in air atmosphere.

Medium Components:
Cholera Toxin from *V. cholerae*,
Catalog Number C8052
DMEM/F12, Catalog Number 51448C
Insulin Solution, Catalog Number I9278
Epidermal Growth Factor, Catalog Number E9644
50 μM Hydrocortisone Solution,
Catalog Number H6909
Horse Serum, Catalog Number H1270

The cryoprotectant medium used is  $1\times$  Cell Freezing Medium-DMSO, Catalog No. C6164.

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

# Biosafety Level: 1

This cell line is not known to harbor an agent known to cause disease in healthy adult humans. Handle as a potentially biohazardous material under at least Biosafety Level 1 containment. The parental cell line, MCF10A, was obtained from ATCC. All animal products used in the preparation of the knockout line and maintenance of both, parental and knockout clone, have been screened negative by 9CFR for adventitious viral agents. Cell lines derived from primate lymphoid tissue may fall under the regulations of 29 CFR 1910.1030 Bloodborne Pathogens. Appropriate safety procedures are recommended to be used when handling all cell lines, especially those derived from human or other primate material. Detailed discussions of laboratory safety procedures have been published. 9-12

## **Preparation Instructions**

Complete Medium: To make the complete growth medium combine the following:

1 L of DMEM/F12 (Catalog Number 51448C) 50 mL of horse serum (Catalog Number H1270) 29 mL of 50  $\mu$ M Hydrocortisone Solution

(Catalog Number H6909)

1.08 mL of Insulin Solution (Catalog Number I9278) 108  $\mu$ L of Cholera toxin solution (1 mg/mL), prepared by dissolving Catalog Number C8052 in sterile water. Store solution at 2–8 °C.

10.8  $\mu$ L of EGF solution (1 mg/mL), prepared by dissolving Catalog Number E9644 in 10 mM acetic acid, followed by 0.2  $\mu$ m filtration. Store the solution in aliquots at –20 °C.

# Storage/Stability

Upon receiving a shipment of frozen cells it is important the end user gives the shipment attention without delay. To ensure the highest level of viability, thaw the vial and initiate the culture as soon as possible upon receipt. If upon arrival, continued storage of the frozen culture is necessary, it should be stored in liquid nitrogen vapor phase and not at –70 °C. Storage at –70 °C will result in loss of viability.

<u>Precaution</u>: It is recommended that protective gloves and clothing always be used, and a full face mask always be worn when handling frozen vials. It is **important to note that some vials leak when submersed in liquid nitrogen** and will slowly fill with liquid nitrogen. Upon thawing, the conversion of the liquid nitrogen back to the gas phase may result in the rapid expansion of the vessel, potentially blowing off its cap with dangerous force creating flying debris.

At the time a cell line is ordered, end users should also consider the culture conditions for the new cell line and make sure the appropriate medium will be available when the cells arrive.

#### **Procedure**

#### Thawing of Frozen Cells

- Thaw the vial by gentle agitation in a 37 °C water bath. To reduce the possibility of contamination, keep the O-ring and cap out of the water. Thawing should be rapid (~2 minutes).
- Remove the vial from the water bath as soon as the contents are thawed, and decontaminate by dipping in or spraying with 70% ethanol. All of the operations from this point on should be carried out under strict aseptic conditions.
- 3. Transfer the vial contents to a centrifuge tube containing 9.0 mL of Complete Medium and spin at  $\sim$ 125  $\times$  g for 5–7 minutes.
- 4. Resuspend cell pellet with the Complete Medium and dispense into a 25 cm² or a 75 cm² culture flask. It is important to avoid excessive alkalinity of the medium during recovery of the cells. It is suggested, prior to the addition of the vial contents, the culture vessel containing the Complete Medium be placed into the incubator for at least 15 minutes to allow the medium to reach its normal pH (7.0–7.6) and temperature (37 °C).
- 5. Incubate the culture at 37 °C in a suitable incubator. A 5% CO<sub>2</sub> in air atmosphere is recommended for the Complete Medium.

#### Subculturing Procedure

Volumes used in this procedure are for a 75 cm<sup>2</sup> flask; proportionally reduce or increase volume of dissociation medium for culture vessels of other sizes.

- 1. Remove and discard culture medium.
- 2. Briefly rinse the cell layer with Trypsin-EDTA solution (Catalog No. T3924).
- Add 2.0–3.0 mL of Trypsin-EDTA solution to the flask and incubate at 37 °C for 15 minutes. This should remove the cells from the cultureware and yield single cells.
- 4. When cells are detached, add 6.0–8.0 ml of Complete Medium and aspirate cells by gentle pipetting.
- Add appropriate aliquots of the cell suspension into new culture vessels.
   Subcultivation Ratio: 1:3 to 1:6
- 6. Incubate cultures at 37 °C.

<u>Note</u>: MCF10A cells require longer time for trypsin digestion than what is typical. More information on enzymatic dissociation and subculturing of cell lines is available in the literature.<sup>13</sup>

#### References

- Zioncheck, T.F. et al., Generation of an active protein-tyrosine kinase from lymphocytes bv proteolysis. J. Biol. Chem., 263, 19195-19202 (1988).
- 2. Coopman, P.J., and Mueller, S.C., The Syk tyrosine kinase: a new negative regulator in tumor growth and progression. Cancer Let., **241**, 159-173 (2006).
- Moroni, M. et al., Progressive loss of Syk and abnormal proliferation in breast cancer cells. Cancer Res., 64, 7346–7354 (2004).
- Sung, Y.M. et al., Tumor Suppressor Function of Syk in Human MCF10A *In Vitro* and Normal Mouse Mammary Epithelium *In Vivo*. PLoS ONE, 4 (10), e7445 (2009). doi: 10.1371/journal.pone.0007445.
- Mahabeleshwar, G.H., and Kundu, G.C., Syk, a protein-tyrosine kinase, suppresses the cell motility and nuclear factor kappa B-mediated secretion of urokinase type plasminogen activator by inhibiting the phosphatidylinositol 3'-kinase activity in breast cancer cells. J Biol Chem., 278, 6209–6221 (2003).
- Larive, R.M. et al., Phosphoproteomic analysis of Syk kinase signaling in human cancer cells reveals its role in cell–cell adhesion Syk kinase signaling in human cancer cells. Oncogene, 28, 2337-2347 (2009).
- 7. Ma, L. et al., The relationship between methylation of the Syk gene in the promoter region and the genesis of lung cancer. Clin. Lab., **56** (9-10), 407-416 (2010).
- Toyama, T. et al., Reduced expression of the Syk gene is correlated with poor prognosis in human breast cancer. Cancer Lett., 189 (1), 97-102 (2003).

- Centers for Disease Control (1999), Biosafety in Microbiological and Biomedical Laboratories Human Health Service Publication No. (CDC) 93-8395. U.S. Dept. of Health and Human Services; 4th Edition U.S. Government Printing Office Washington D.C. The entire text is available online at <a href="https://www.cdc.gov/od/ohs/biosfty/bmbl4/bmbl4toc.htm">www.cdc.gov/od/ohs/biosfty/bmbl4/bmbl4toc.htm</a>
- Fleming, D.O. et al., (1995) Laboratory Safety: Principles and Practice. Second edition, ASM press, Washington, DC.
- Hay, R.J. et al., eds. (1992), ATCC Quality Control Methods for Cell Lines. 2<sup>nd</sup> edition, Published by ATCC.
- 12. Caputo, J.L., Biosafety procedures in cell culture. J. Tissue Culture Methods, **11**, 223-227 (1988).
- 13. Freshney, R.I., Chapter 10 in Culture of Animal Cells, a manual of Basic Technique by, 3rd edition, published by Alan R. Liss, (NY, NY: 1994).

Additional product and technical information can be obtained from the catalog references and the Sigma Life Science Website (www.wherebiobegins.com/biocells).

CompoZr is a registered trademark of Sigma-Aldrich® Biotechnology LP and Sigma-Aldrich Co. ATCC is a registered trademark of American Type Culture Collection.

CRL-10317 is a trademark of American Type Culture

CRL-10317 is a trademark of American Type Culture Collection.

Please see the enclosed Label License Agreement (LLA) for further details regarding the use of this product. The LLA is also available on our website at www.wherebiobegins.com/biocells

These cells are distributed for research purposes only. Sigma Life Science requires that individuals contemplating commercial use of any cell line first contact us to negotiate an agreement. Third party distribution of this cell line is prohibited.

GW,ADM,MAM 06/11-1