3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

LYN A, active, GST-tagged, human PRECISIO[®] Kinase recombinant, expressed in *Sf*9 cells

Catalog Number **SRP5277** Storage Temperature –70 °C

Synonyms: JTK8, FLJ26625

Product Description

LYN A is a member of the Src-family of protein kinases that participate in signaling pathways of a variety of cell surface receptors and are localized to the cytoplasmic side of the plasma membrane. LYN A is biosynthetically transported to the plasma membrane via the Golgi pool of caveolin along the secretory pathway. The trafficking of LYN A from the Golgi apparatus to the plasma membrane is inhibited by deletion of the kinase domain but not by kinase inactivation. The kinase domain of LYN A plays a role in LYN A trafficking besides catalysis of substrate phosphorylation. cDNA microarray analysis revealed that LYN A may be involved in the changes that are related to tumor pathogenesis.²

Recombinant full-length human LYN A was expressed by baculovirus in *Sf*9 insect cells using an N-terminal GST-tag. The gene accession number is NM_002350. It is supplied in 50 mM Tris-HCl, pH 7.5, 150 mM NaCl, 10 mM glutathione, 0.1 mM EDTA, 0.25 mM DTT, 0.1 mM PMSF, and 25% glycerol.

Molecular mass: ~81 kDa

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The product ships on dry ice and storage at -70 °C is recommended. After opening, aliquot into smaller quantities and store at -70 °C. Avoid repeated handling and multiple freeze/thaw cycles.

Figure 1.

SDS-PAGE Gel of Typical Lot:

≥70% (SDS-PAGE, densitometry)





Figure 2.
Specific Activity of Typical Lot: 434–657 nmole/min/mg



Procedure

Preparation Instructions

Kinase Assay Buffer – 25 mM MOPS, pH 7. 2, 12.5 mM glycerol 2-phosphate, 20 mM MgC1₂, 12.5 mM MnC1₂, 5 mM EGTA, and 2 mM EDTA. Just prior to use, add DTT to a final concentration of 0.25 mM.

Kinase Dilution Buffer – Dilute the Kinase Assay Buffer 5-fold with a 50 ng/μL BSA solution.

Kinase Solution – Dilute the active LYN A ($0.1 \mu g/\mu L$) with Kinase Dilution Buffer to the desired concentration. Note: The specific activity plot may be used as a guideline (see Figure 2). It is recommended the researcher perform a serial dilution of active LYN A kinase for optimal results.

10 mM ATP Stock Solution – Dissolve 55 mg of ATP in 10 mL of Kinase Assay Buffer. Store in 200 μ L aliquots at –20 °C.

 $\gamma\text{-}^{33}\text{P-ATP}$ Assay Cocktail (250 $\mu\text{M})$ – Combine 5.75 mL of Kinase Assay Buffer, 150 μL of 10 mM ATP Stock Solution, 100 μL of $\gamma\text{-}^{33}\text{P-ATP}$ (1 mCi/100 μL). Store in 1 mL aliquots at –20 °C.

Substrate Solution – Dissolve the protein substrate in distilled water at a final concentration of 1 mg/mL.

1% phosphoric acid solution – Dilute 10 mL of concentrated phosphoric acid to a final volume of 1 L with water.

Kinase Assay

This assay involves the use of the ³³P radioisotope. All institutional guidelines regarding the use of radioisotopes should be followed.

- 1. Thaw the active LYN A, Kinase Assay Buffer, Substrate Solution, and Kinase Dilution Buffer on ice. The γ -33P-ATP Assay Cocktail may be thawed at room temperature.
- 2. In a pre-cooled microcentrifuge tube, add the following solutions to a volume of 20 μ L:

10 μL of Kinase Solution

 $5 \,\mu\text{L}$ of Substrate Solution

5 μL of cold water (4 °C)

- 3. Set up a blank control as outlined in step 2, substituting 5 μ L of cold water (4 °C) for the Substrate Solution.
- 4. Initiate each reaction with the addition of 5 μ L of the γ - 33 P-ATP Assay Cocktail, bringing the final reaction volume to 25 μ L. Incubate the mixture in a water bath at 30 °C for 15 minutes.
- After the 15 minute incubation, stop the reaction by spotting 20 μL of the reaction mixture onto an individually precut strip of phosphocellulose P81 paper.

- Air dry the precut P81 strip and sequentially wash in the 1% phosphoric acid solution with constant gentle stirring. It is recommended the strips be washed a total of 3 times of ~10 minutes each.
- 7. Set up a radioactive control to measure the total γ - 33 P-ATP counts introduced into the reaction. Spot 5 μ L of the γ - 33 P-ATP Assay Cocktail on a precut P81 strip. Dry the sample for 2 minutes and read the counts. Do not wash this sample.
- 8. Count the radioactivity on the P81 paper in the presence of scintillation fluid in a scintillation counter.
- 9. Determine the corrected cpm by subtracting the blank control value (see step 3) from each sample and calculate the kinase specific activity.

Calculations:

1. Specific Radioactivity (SR) of ATP (cpm/nmole)

SR =
$$\frac{\text{cpm of 5} \mu \text{L of } \gamma^{-33} \text{P-ATP Assay Cocktail}}{\text{nmole of ATP}}$$

cpm – value from control (step 7) nmole – 1.25 nmole (5 μ L of 250 μ M ATP Assay Cocktail)

2. Specific Kinase Activity (SA) (nmole/min/mg)

nmole/min/mg =
$$\Delta cpm \times (25/20)$$

SR \times E \times T

SR = specific radioactivity of the ATP (cpm/nmole ATP) ∆cpm = cpm of the sample – cpm of the blank (step 3) 25 = total reaction volume

20 = spot volume

T = reaction time (minutes)

E = amount of enzyme (mg)

References

- Kasahara, K. et al., Trafficking of Lyn through the Golgi caveolin involves the charged residues on alphaE and alphal helices in the kinase domain. J. Cell Biol., 165(5), 641-52 (2004).
- Baran, C.P. et al., The inositol 5'-phosphatase SHIP-1 and the Src kinase Lyn negatively regulate macrophage colony stimulating factor-induced Akt activity. J. Biol. Chem., 278(40), 38628-36 (2003).

PRECISIO is a registered trademark of Sigma-Aldrich Co. LLC.

RC.MAM 12/12-1