

ProductInformation

CALF THYMUS DNA, Ultrapure

Product No. **D 4764**

Lot 124K6103

Storage: 2-8 °C

Product Summary

Approximate Average Size: >13 kb

A_{260}/A_{280} : 1.9

Unit Definition

One unit will yield an A_{260} of 1.0 in 1.0 ml of 1 mM NaCl, 1 mM Tris-HCl pH 7.5, and 1 mM EDTA (1cm light path).
One mg of DNA is equivalent to approx. 20 units.

Description

This DNA has been purified by equilibrium buoyant density gradient ultracentrifugation in cesium chloride. The DNA is collected from the gradient and dialyzed against 1 mM NaCl, 1 mM Tris-HCl, pH 7.5, and 1 mM EDTA. The solution was then adjusted to a concentration of 8.93 A_{260} units of DNA per ml and lyophilized.

Size Determination by Agarose Gel Electrophoresis

Agarose gel electrophoresis was performed on a 0.7% agarose gel using 1X TAE buffer. The gel was electrophoresed for 3 hours at 70 volts. The reference DNA used for size determination was single cut λ DNA (D 4917) and this lot was determined to have an average size of >13 Kb.

References:

1. Burgi, E., and Hershey, A., *Biophys. J.*, **3**, 309 (1963).
2. Mandel, M. and Marmur, J., *Methods Enz.*, Vol. XII, Part B, 195 (1968).

4/05

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.